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Hear now. And always.

Notice of update to MRI Guidelines and Patient Information

Cochlear is providing notice that updates to the the following documents are in process but were not available at the time of this shipment:

- Document Number D774756: Cochlear™ Nucleus® Implants Magnetic Resonance Imaging (MRI) Guidelines
- Document Number 563636: Nucleus CI532 cochlear implant Patient Information
- Document Number 551440: Nucleus CI522 cochlear implant Patient Information
- Document Number 491848: Nucleus CI512 cochlear implant Patient Information
- Document Number 216687: Nucleus CI24RE (CA) cochlear implant Patient Information
- Document Number 216687: Nucleus CI24RE (ST) cochlear implant Patient Information
- Document Number 264715: Nucleus CI422 cochlear implant Patient Information
- Document Number 415754: Nucleus CI24REH cochlear implant Patient Information

Effective January 6, 2020, the Cochlear Nucleus Implant Bandage and Splint Kit for MRI (MRI Kit) is no longer offered by Cochlear in the United States.

The MRI Kit was offered for use with certain implant types shown below to provide access to MR scan at 1.5 Tesla. Based on these changes, recipients with these implants can obtain MR scans with the following considerations:

- The implant magnet must be removed for 1.5 T and 3 T scans.
- Follow all other conditions for safe scanning given in the MRI Safety Information in the MRI Guidelines documentation.

Consideration and steps for magnet removal are found in the implant Physician Guides and MRI Guidelines.

What implant models are affected by this change?

Profile Series: CI512, CI522, CI532

Freedom Series: CI422, CI24REH, CI24RE (CA), CI24RE (ST)

For any further questions about the changes within the MRI Guidelines, please call MRI Information Line at 1-866-210-9217 or contact your local Cochlear representative.

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Tips for HiRes 90K Magnet Removal—Replacement Procedure

Disclaimer

Advanced Bionics® does not endorse or recommend any procedures other than those stated in the HiRes 90K® Surgeon's Manual for HiFocus® Helix and HiFocus® 1j Electrodes, 9055112-002 and the HiRes 90K® Surgical Video, CI-8167.

We are grateful to

Ernesto Diaz-Ordaz, MD
Buffalo Hearing & Speech Center
Sisters of Charity Hospital
Buffalo, New York

The surgical tips discussed below are available in a PowerPoint presentation from Boston Scientific.

Things to consider when planning a magnet removal—replacement procedure:

- Use bipolar electrocautery only.
- Have a plastic HiRes 90K Mock-up, CI-4425, to help determine the position of the implanted device and assist in marking the proposed incision line.
- Incision is posterior to the implant coil-antenna. Avoid placing the incision directly over the implanted device.
- When working with the implant magnet or replacement magnet, CI-1412, consider using non-magnetic tools and/or alternative techniques (for example, titanium retractors or sutures to retract the skin flap).
- Have a non-metallic gimmick available to help position the replacement magnet (for example, the plastic "filler" from a suction line).
- Use a headpiece, placed in a sterile glove, to test for proper magnet orientation ("+" side up) and adequate headpiece retention.
- Close the wound and cover with a light-pressure, mastoid dressing.
- Have a back-up HiRes 90K device available in case the implant is damaged during surgery.

Magnet-Removal Procedure

Perform a minimal shave posterior to the device magnet-coil.

Use a plastic HiRes 90K Mock-up, CI-4425, to help determine the position of the implanted device and assist in marking the proposed incision line.

Incision


Use bipolar electrocautery only.

Raise the skin-muscle flap over the internal magnet-coil.


Use a freer elevator to access and lift the original magnet (which should be attracted to the freer elevator) from the implant.

Since most surgical tools are metallic, manipulation of the replacement magnet and its insertion into the implant can be problematic.


Use sutures to retract the skin flap.

Manually place the replacement magnet, CI-1412, (“+” side up) in the recess of the implant.

A non-metallic gimmick—the plastic “filler” from a suction line—was used to seat the magnet in the recess under the silicone lip.

Replacement Magnet, CI-1412, in place.

Use a headpiece, placed in a sterile glove, to test for proper magnet orientation (“+” side up) and adequate headpiece retention.

Close the wound and cover with a light-pressure, mastoid dressing.

ADVANCED BIONICS GmbH | Max-Eyth-Str. 20 | DE-70736 Fellbach

Ansprechpartner: Deike Markgraf
Telefon: +49.172.184.06.29
Fax: +49.711.510.70.571
Web: advancedbionics.com
E-Mail: Deike.Markgraf@AdvancedBionics.com

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A. Schwaiblmair

Advanced Bionics GmbH, Max-Eyth-Str. 20, 70736 Fellbach-Oeffingen, Germany
+49.711.510.70.570 – info.de@advancedbionics.com

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Cochlear implant magnet dislocation after MRI: surgical management and outcome

Martin Leinung¹ · Andreas Loth¹ · Maximilian Gröger¹ · Iris Burck² · Thomas Vogl² · Timo Stöver¹ · Silke Helbig¹

Received: 4 May 2019 / Accepted: 22 January 2020
© Springer-Verlag GmbH Germany, part of Springer Nature 2020

Abstract

Purpose An increasing number of cochlear implant (CI) users is examined by magnetic resonance imaging which may cause the displacement of the implant magnet. This complication prevents the usage of the external processor and has to be treated surgically in most cases. The purpose of this study is to analyze the results of the surgical intervention and the consequences for the CI recipients.

Methods The retrospective study was conducted at a tertiary referral center. From the patient care records between October 2014 and July 2018, 9 cases were reviewed that had undergone MRI after cochlear implantation and had experienced magnet displacement.

Results Nine patients from 9 to 74 years of age were identified with MRI-induced magnet displacement. Implants of different manufacturers were affected (8 × Cochlear[®], 1 Advanced Bionics[®]) but did not include the latest 3 T MR conditional product generation. The patients reported pain, swelling, redness above the implant and/or a noticeably dislocated magnet. One-third of the MRI examination were conducted in external radiological sites without any precautions such as a compression bandage. Surgical magnet repositioning was successful in all but one case with postoperative implant infection and consecutive explantation. In total, the patient was unable to use his CI for 420 days (1.2 years) after the MRI examination. The remaining eight patients averaged 29 days between MRI-related magnet dislocation and CI re-activation.

Conclusions The present study shows that in the majority of cases a surgical magnet reposition is possible without complications, and thus the time of nonuse of the CI is usually low. Nevertheless, there is a risk that in individual cases significant medical, functional, social and economic consequences for patients may occur. The presented data demonstrate that the indication to perform MRI scans in CI users needs to be further critically considered. An attentive, critical assessment of an MRI indication by both the initiating physician (usually not an ENT specialist) and the performing radiologist is mandatory.

Keywords Cochlear implant · Magnet displacement · Magnetic resonance imaging

Introduction

The cochlear implant (CI) has been established as the gold standard in the functional rehabilitation of severe hearing loss and deafness over the last decades [1]. Originally conceptualized as a measure to support the lip-reading of affected patients, meanwhile, the CI treatment aims for open

speech understanding, verbal communication at the phone and good understanding in noise [2, 3]. In addition, a number of indication changes have taken place in recent years. These include amongst others double-sided (bilateral) CI surgery [4, 5], CI treatment for unilateral deafness [6], and especially CI therapy for elderly and very elderly (geriatric) patients [7, 8]. Thus, the CI is the success story in neurobiotics. These developments have led to a steady increase in the number of patients treated with a CI over the last decades. Worldwide, it is estimated that about 500,000 subjects were supplied with a CI [9].

Parallel to this, there has been a significant increase on radiologic imaging in the last decades in the context of differential diagnostic examinations, in particular with the help of MRI technology [10]. Numerous indications

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¹ Department of Otorhinolaryngology, Head and Neck Surgery, University Hospital Frankfurt, Theodor-Stern-Kai 7, 60590 Frankfurt am Main, Germany

² Institute of Diagnostic and Interventional Radiology, University Hospital, Frankfurt, Germany

have been established here in addition to the examination of the neurocranium. These include the examination of the neural structures of the spine (cervical, thoracic, lumbar spine) or of joint structures of the musculoskeletal system (knee, hand, shoulder, etc.). In an inpatient setting the number of MRIs per year in Germany more than doubled from ca 1,000,000 in 2005 to 2,034,000 in 2017 [11]. Considering the effects of the aging population and extended CI indications, there is a significant need for performing MRI scans even in patients who have already received a CI.

While in the early days following the introduction of CI care, there was a clear contraindication regarding MR imaging in CI patients, this situation has changed in recent years. Meanwhile, all four major manufacturers (alphabetically: Advanced Bionics, Cochlear Ltd., MED-EL, Oticon Medical) offer current CI products which are approved for MRI or "MRI-conditional".

In addition to an inevitable artifact due to the metallic components of the CI, MRI scans of a cochlear implant endanger the demagnetization or polarity reversal of the magnet, pain and displacement of the magnet out of its housing [12–18]. Magnet displacement can lead to problems with fixing the speech processor, discomfort or even skin necrosis and infection and needs to be treated by surgical magnet reposition or exchange [19, 20]. Three manufacturers now use rotatable magnet components to minimize the risk of demagnetization and displacement [17, 21]. First reports in a 3 T MRI were promising as the scan had been performed without any pain sensation or other adverse effects [22]. Generally, manufacturers recommend a tight wrap bandage over the implant to prevent magnet dislocation during MRI [23, 24]. The practical implementation of this pressure bandage, however, often presents a challenge due to individual head size and shape, varying amount of hair and soft tissue, so that despite this protective measure a magnet dislocation can occur during the MRI examination [25–28].

In case of a dislocation a surgical repositioning is required in almost every case. The following steps after the suspected diagnosis of magnet displacement are radiological verification, surgical intervention, postoperative wound healing and re-activation of the implant. In the meantime, the CI processor unit usually cannot be worn. So far, there are only a few studies available on the time of non-usage of the CI due to magnet dislocation and possible complications due to MR-related magnet dislocation. The aim of the study presented here was, therefore, to examine the results of the surgical procedure, the resulting periods of non-use of the CI ("offline time") and the collection of possible complications in connection with the MR-induced magnet dislocation of a CI.

Patients and methods

For this retrospective study consent of the local ethics committee of the medical faculty of our institution was obtained (No. 391/17). The study took the collection of data from patient care records between October 2014 and July 2018 into account. From the available patient records, cases were identified that had undergone MRI after cochlear implantation and had experienced magnet displacement. The primary aim of the study was to record the results of the surgical intervention and the consequences for the patient. This is in particular the period of inoperability of the implant ("offline time") with the consequent result of transient functional deafness. For this purpose, the following data were collected: Patient's age in years at the time of the MRI scan, cochlear implant (manufacturer and type), period between CI implantation and MRI performed (in months), indication for MRI examination (organ), MRI field strength in Tesla, use of a head bandage during MRI, site of MRI examination (internal/external), symptoms after MRI as an indication of magnet dislocation (pain, swelling, redness or palpable displacement, inability to wear the external processor), method for confirming magnet dislocation (conventional X-ray, CT, cone beam CT); time between MRI and magnet repositioning (in days), time period between magnet repositioning and CI activation (in days), presence of wound infection due to magnet repositioning, total period of nonuse of CI ("offline time") in days.

Results

A total of nine patients were identified who had undergone proven magnet dislocation during an MRI scan. The result overview is shown in Table 1. Among the affected patients were 7 women and 2 men. The age of the affected patients at the time of the MRI examination averaged 34 years (median, min 9, max 74 years).

The time from CI surgery to MRI varied greatly between 0.3 and 108 months (median 23 months). Of the 9 dislocation cases, 4 were lumbar spine examinations for orthopedic symptoms and 2 cranial MRIs for neurological symptoms. In one case each knee, wrist and thoracic spine were scanned. Five out of nine MRI examinations with subsequent dislocation were performed externally, outside of our facility, while four were performed internally. All internally scanned patients received a compression bandage including a stiffening element (plastic card or plastic cylinder), which was applied by an ENT specialist prior to the scan. Of the five external patients, two individuals

Table 1 Summary of clinical patient data

Patient no	1	2	3	4	5	6	7	8	9	Median
Implant										
Manufacturer	CL	CL	CL	CL	CL	CL	AB	CL	CL	
Type	CI 512	CI 24RE	CI 512	CI 512	CI 24RE	CI 532	HiRes 90 k	CI 512	CI 512	
MRI scan										
Age at MRI [years]	37	9	34	72	34	74	21	11	54	34
Time CI to MRI [months]	51	88	34	10	108	8	0.3	23	2	23
Anatomic region scanned	lumbar	cranial	cranial	lumbar	hand	lumbar	lumbar	thorax	knee	
Magnet field strength [T]	1.5	uk	uk	1.5	uk	uk	1.5	1.5	1.5	
Bandage during MRI [yes●/no○]	●	uk	●	●	○	uk	●	●	●	
MRI site [i: intern/e: extern]	i	e	e	i	e	e	i	i	e	
Magnet dislocation										
Complaints [yes●/no○]										
Pain	○	●	●	○	●	○	○	○	●	
Swelling	●	○	●	●	○	●	○	○	●	
Redness	○	●	○	○	○	○	○	○	●	
Palpable displacement	○	●	○	○	●	○	●	○	○	
Inability to wear the CI processor	●	●	●	○	●	●	●	○	●	
Confirming imaging modality	CT	X-ray	CBCT	X-ray	CT	CBCT	CBCT	X-ray	X-ray	
Time to magnet replacement [days]	8	67	7	7	4	8	5	6	9	7
Time to CI reactivation [days]	11	28	17	31	7	412	27	20	20	20
Total CI „offline“ time [days]	19	95	24	38	11	420	32	26	29	29
Infection [yes●/no○]	○	○	○	○	○	●	○	○	○	

CL Cochlear Limited, AB Advanced Bionics, uk unknown

were protected with a bandage, one external patient did not receive any precautionary measure against magnet dislocation and two patients could not recall it. In five out of nine cases a field strength of 1.5 T was used to perform the MRI examination. In four cases the applied field strength is not known due to the external execution of the examination. During or after the study, affected patients reported single or combined symptoms that suggested magnet dislocation. The symptoms included pain, swelling, redness above the implant and/or a dislocated magnet felt by the patient. Magnet dislocation was confirmed in four out of nine cases by conventional radiographs, twice by CT and three times via CBCT (see Fig. 1).

The survey of the cochlear implants used showed that in eight cases an implant from Cochlear Ltd. (5 × CI512, 2 × CI24 RE, 1 × CI532) and in one case an implant from Advanced Bionics (HiRes90K) was affected by magnet displacement.

After radiographic verification surgical reposition of the magnet was recommended to the affected patients. All but one patient made immediate use of this offer so that the intervention to reposition took place in eight of nine patients within 4–9 days (median 7 days). Due to the individual patient's request, the surgery was done in one case only after 67 days (patient no. 2).

In all nine cases the surgical magnet repositioning was performed successfully under general anesthesia and antibiotic prophylaxis. The postoperative wound healing was uneventful in eight out of nine cases. In one case (patient no. 6) a chronic wound healing disorder occurred with bacterial superinfection and repeated systemic antibiotic treatments. Finally, the explantation of the device was inevitable. The re-implantation could only take place after a further 3 months, so that the patient was able to actively use the cochlear implant only 420 days or 1.2 years after MRI-induced magnet dislocation. The remaining eight patients averaged 29 days between MRI-related magnet dislocation and CI re-activation (min: 11 days, max: 95 days).

Discussion

The major achievements of hearing rehabilitation using a cochlear implant have led to a significant dissemination of the method. However, it also brought up new challenges as its impact on the diagnosis and therapy of hearing-independent diseases has to be considered: The application of MRI diagnostics necessitates particular demands on a medical device with respect to patient safety [21]. This is especially


Fig. 1 Modified Stenvers view radiography before (a) and after (b) magnet displacement. The implant is a Nucleus® CI24RE (1) with a contour advance electrode. Before displacement, the magnet (2) is located in the geometrical centre of the antenna coil (3). In panel b the magnet is clearly shifted outwards in relation to the antenna coil

(white arrow). The other implant components are the reference electrode (4), the mastoidal sigmoid course of the electrode lead (5) and the intracochlear electrode array (6). Between the images the ventilation tube had been removed (black arrow) and a contralateral CI had been implanted (astersik)

true for a CI, because in addition to the pure consideration of a potential damage to the implant or even the patient, these groups of patients additionally suffer from the nonuse of a CI by the temporary loss of hearing and communication. To our best knowledge dislocation, this is the first study with focus on the "offline time", i.e. the duration of the patient impairment resulting from magnet dislocation after MR imaging.

Precautions against magnet dislocation

Table 2 shows a summary of the mandatory and optional measures for MRI imaging of CI users for the field strengths of 1.5 and 3 T commonly used today as

recommended by the different manufacturers [23, 29–33]. Very early implant types were not included for reasons of clarity; the selection of implants made covers > 95% of the CI recipients implanted at the clinic site since 1987. The implant generations of the company Advanced Bionics illustrate the technological progress: While older implants were not even approved for 3 T MRI, the penultimate family of devices enabled 3 T imaging, at least under safety measures (head bandage after magnet removal). And finally, the latest product allows both 1.5 T and 3 T imaging completely without any safety measures, provided that general conditions regarding the maximum specific absorption rate (SAR) and the position of the head are

Table 2 Requirements for 1.5 T and 3 T MR imaging for the most common CI types currently in use

manufacturer and device	1,5T MRI		3T MRI		● events of magnet dislocation	cumulative incidence of magnet dislocation
	bandage mandatory	magnet explanation	bandage mandatory	magnet explanation		
AB HiRes90K (Advantage)® HiRes Ultra® HiRes Ultra 3D®	mandatory	optional	mandatory	mandatory	●	dislocation cases: 1 implants till 2018: 190 incidence: 0.53%
	mandatory	optional	mandatory	mandatory	●●●●●	
	optional	optional	optional	optional	●	
Cochlear Ltd. Nucleus Series CI24RE® CI512® Nucleus Profile Series CI522® CI532® Nucleus Profile Plus Series CI612® CI622® CI632®	mandatory	optional	mandatory	mandatory	●●●●●	dislocation cases: 8 implants till 2018: 1,075 incidence: 0.74%
	optional	optional	optional	optional	●	
	mandatory	magnet not removable	mandatory	mandatory		
	optional	optional	optional	optional		
MED-EL C40®, C40+® PULSARci100® SONATATI100® Mi1000 CONCERTO (PIN)® Mi1200 SYNCHRONY (PIN/ST)® Mi1210 Synchrony ST® Mi1250 SYNCHRONY 2 (PIN)®	mandatory	magnet not removable	mandatory	mandatory		dislocation cases: 0 implants till 2018: 1,299 incidence: 0.00%
	optional	optional	optional	optional		
	mandatory	optional	mandatory	mandatory		
	mandatory	optional	mandatory	mandatory		

AB Advanced Bionics, OM Oticon Medical

Red dots symbolize magnet dislocation events in our study population

observed. A compression bandage may be an option if wound healing, infection or trauma make destabilization of the implant site appear possible. As a rule of thumb, a wrapping bandage is recommended until 6 months after implantation [33]. Depending on the underlying medical issue, elective removal of the magnet may be useful to limit implant-related artifacts [19]. Red dots in Table 2 represent individual cases of dislocation in our study population. Here we see a quantitative predominance of eight dislocations for implants from the manufacturer Cochlear, followed by one dislocation for Advanced Bionics implants, whereas we did not find a single dislocation for a MED-EL implant. This pattern was also evident in the patient groups presented by other authors [19]. The number of adverse events is generally determined by the probability of occurrence of the event and the number of procedures performed. In our clinic, cochlear products are implanted almost six times more frequently than AB implants. Assuming that the likelihood of MRI examination is the same for each type of implant, the probability of occurrence of magnet dislocation can be calculated as shown in the right column of Table 2. The resulting cumulative incidences of 0.53% for Advanced Bionics and 0.74% for Cochlear are comparably low.

MED-EL's flawless result with not a single magnet dislocation has two reasons: In previous product generations up to the Mi1000 CONCERTO® the magnets were fixed. This had the disadvantage that the image artifact could not be reduced by removing the magnets, but effectively prevented dislocation. The second reason is that MED-EL introduced a sufficient MRI conditional technology with the first SYNCHRONY® implant in 2014, several years ahead of the other manufacturers. Meanwhile, all three of the most widely used manufacturers have 3 T-compatible implants in their portfolio, which promise a comfortable and safe imaging procedure. To our best knowledge, the latest generations of implants have not yet experienced a dislocation event. In the long term, therefore, the problem of magnet dislocation caused by MRI should fade into the background as these implants become more widespread.

Until then, however, radiologists and ENT physicians will have to deal with the problem that in Germany alone a large part of the approximately 100,000 CI users have been fitted with implants that are only partially MRI-compatible and will increasingly undergo MRI examinations of any part of the body [11, 23]. The age of our patient group was between 9 and 74 years (median 34 years). Considering the period between the initial cochlear implantation and the MRI examination, a substantial variance is noticeable: The minimum period was 0.3 months (11 days, patient No. 2) up to 108 months (9 years, patient No. 5). These data show that a dislocation during MR imaging can occur at any age and any time after CI surgery.

A compression bandage is the most often recommended countermeasure against magnet displacement. Cochlear® recommends a "splint" directly above the internal magnet in addition to the compression bandage since 2011 [34]. Force measurements conducted by our research group have shown, that a stiffening element increases the stabilization of the magnet significantly compared to a compression bandage alone (publication in preparation). A compression bandage was applied in six of nine patients and at least the four internally scanned patients received an additional stiffening element. But despite all the precautions carefully taken the magnets were displaced during the MRI scan. Several authors had the same unfortunate experience [18, 20, 34]. The benefit of a compression bandage is difficult to quantify: Young et al. experienced one displacement event in 18 MRI scans. This would result in a failure rate of the CI bandage of 5.6% [18], but the number of patients is small. Alone four different implant types were included in the study. Therefore, it is impossible to break down the numbers for each implant. Other factors such as the position of the patient in the magnetic field of the MR scanner and the scanning protocol itself will most likely have an additional impact. In addition, a presumably relevant number of unreported cases of externally performed MRI examinations, which are not known to the implantation clinic, can be assumed.

The retrospective questionnaire study of Grupe et al. [35] revealed a prevalence of 9.8% that a CI patient had one or more MRI scans in the past. 20 patients (61%) had received a compression bandage during MR imaging (which is comparable to the 6 of 9 = 67% of the patients in our study). 1 of 3 magnet dislocation events occurred without a bandage. The resulting bandage failure rate of $1/20 = 5\%$ is similar to the 5.6% of Young et al. (see above) [18]. More reliable numbers on the incidence of MR induced magnet displacement and the benefit of compression bandages could be obtained if many more asymptomatic patients were interviewed, e.g. through the collaboration of implantation clinics. A national implant register could make a meaningful contribution to this, provided that the database is designed in such a way as to enable the long-term assignment of complications to individual pseudonymized patients and their implants.

A possible approach for improvement could be to carry out MRI examinations of CI patients primarily at specialized centers (most likely implantation clinics). Although this by no means excludes magnetic dislocations (as the number of four internal cases in our patient group shows) such a model offers several advantages: ENT specialists on site can clarify more quickly which precautionary measures are really necessary in individual cases. Trained radiologists are familiar with the scanning protocols to be followed and are alert to the clinical signs of magnetic dislocation. And finally, a higher number of cases is often a successful way to minimize complication rates.

Diagnosis

The first challenges after a magnet dislocation are the early suspected diagnosis and subsequent confirmation of the event. Our patients experienced pain, redness, swelling over the implant or a palpable elevation in the area of the dislocated magnet. These symptoms occurred solely or in any possible combination, either during MR imaging or directly following the examination. Most of the patients were not able to wear the external processor anymore and therefore suffered from hearing loss.

From the variability of the reported symptoms, it becomes clear that there are often no definite or clear symptoms of magnet dislocation. One individual (patient no. 8) even did not recall any of the four symptoms and was diagnosed accidentally during a regular checkup in our clinic. One may assume that there is a dark figure of undetected, asymptomatic cases of magnet dislocation. None of the MRI examinations of our nine patients had been indicated by an ENT physician. Since it cannot be assumed that an ENT physician is involved in every MRI imaging of CI patients, the examining radiologist must be particularly sensitized to this topic. The critical and active request of symptoms immediately after the MRI examination, e.g. by means of a standardized questionnaire, could contribute to the immediate assurance of the diagnosis. This would be beneficial as in all cases a radiological method confirmed the suspicion. X-ray, CBCT or CT seem to be equally suitable to detect magnet dislocation, wherefore X-ray is preferred due to minimum radiation exposure. Hassepas et al. recommend to carry out an X-ray in Stenvers projection directly after every MRI in CI users [34].

Surgery and CI offline time

In case of a proven dislocation, the therapy almost always consists of surgical repositioning. An exception is a grade I displacement according to Gubbels and McMenomey (canting of the magnet within its housing, the silicone lip however, is still on the superficial and lateral side of the magnet [25]). In that case a manual reposition from outside is feasible [36]. Every higher grade of displacement


requires a surgical intervention. A common incision is a semicircular cut around the antenna area of the implant as shown in Fig. 2a [34, 37]. Ray et al. found that larger incisions correlate with wound and flap problems [38]. For this reason, Leong et al. described an endoscopic repositioning technique using a smaller access route (B) [20]. We prefer an open surgical approach in combination with a small semicircular skin incision between the magnet and the antenna above the implant (C). In most cases the procedure was conducted under antibiotic prophylaxis in general anesthesia but could also be performed in analgesedation [37] or local anesthesia [19].

The time between magnet dislocation and reposition surgery was 4–9 days (median 7 days) with the exception of patient No. 2 who asked to delay the surgery for 67 days. But even if this individual is omitted the median stays at 7 days (a median is less sensitive to outliers in smaller sample numbers than the average). Reactivation of the implant took place 7 days after surgery (median). This total offline time of 29 days (median, without patients 2 and 6: 26 days, see below) can be shortened by up to a week if the repositioning is performed under local anesthesia on the same day as the MRI scan as the treatment concept of Hassepas et al. suggests [34].

However, the essential point is the postoperative wound infection of patient No. 6: Despite all local and systemic antibiotic and operative therapy attempts the explantation of the CI was inevitable. Following a 3 months healing period, the patient was subsequently successfully re-treated with a CI and achieved a similar level of speech intelligibility as before the event. Overall, this resulted in a phase of 60 weeks in which the patient could not use his CI (offline time).

This severe complication had been previously reported by other authors as well [19, 20] and in the long run it can be solved. But in the meantime, i.e. for 1.2 years, the patient was significantly impaired and considerably affected in the everyday situation and quality of life. The case undoubtedly represents a worst-case scenario that could only be surpassed if an infant in the middle of its speech development had been affected. However, it illustrates the need to make physicians of all disciplines aware that MRIs for CI patients are still not

Fig. 2 Various skin incision lines (red) for magnet replacement

a matter of course and are associated with relevant risks. An attentive, critical assessment of an MRI indication by both the initiating physician (usually not an ENT specialist) and the performing radiologist is mandatory.

Compliance with ethical standards

Conflict of interest There are no conflicts of interest to disclose by the authors.

Human and animal rights The submitted manuscript does not report data derived from experimental or clinical observations in human or animal subjects.

Informed consent For the same reason informed consent was neither necessary nor applicable.

References

- von Ilberg CA, Baumann U, Kiefer J, Tillein J, Adunka OF (2011) Electric-acoustic stimulation of the auditory system: a review of the first decade. *Audiol Neurotol* 16(Suppl 2):1–30. <https://doi.org/10.1159/000327765>
- Holden LK, Finley CC, Firszt JB, Holden TA, Brenner C, Potts LG, Gotter BD, Vanderhoof SS, Mispagel K, Heydebrand G, Skinner MW (2013) Factors affecting open-set word recognition in adults with cochlear implants. *Ear Hear* 34(3):342–360. <https://doi.org/10.1097/AUD.0b013e3182741aa7>
- Mudery JA, Francis R, McCrary H, Jacob A (2017) Older individuals meeting medicare cochlear implant candidacy criteria in noise but not in quiet: are these patients improved by surgery? *Otol Neurotol* 38(2):187–191. <https://doi.org/10.1097/MAO.0000000000001271>
- Kraaijenka JVC, Ramakers GJJ, Smulders YE, van Zon A, Stegeman I, Smit AL, Stokroos RJ, Hendrice N, Free RH, Maat B, Frijs JHM, Briare JJ, Mylanus EAM, Huinck WJ, Van Zanten GA, Grolman W (2017) Objective and subjective measures of simultaneous vs sequential bilateral cochlear implants in adults: a randomized clinical trial. *JAMA Otolaryngol Head Neck Surg* 143(9):881–890. <https://doi.org/10.1001/jamaoto.2017.0745>
- Goupell MJ, Kan A, Litovsky RY (2016) Spatial attention in bilateral cochlear-implant users. *J Acoust Soc Am* 140(3):1652
- Dorman MF, Natale SC, Butts AM, Zeitler DM, Carlson ML (2017) The sound quality of Cochlear implants: studies with single-sided deaf patients. *Otol Neurotol* 38(8):e268–e273. <https://doi.org/10.1097/MAO.0000000000001449>
- Mosnier I, Bebear JP, Marx M, Fraysse B, Truy E, Lina-Granade G, Mondain M, Sterkers-Artières F, Bordure P, Robier A, Godey B, Meyer B, Frachet B, Poncet-Wallet C, Bouccara D, Sterkers O (2015) Improvement of cognitive function after cochlear implantation in elderly patients. *JAMA Otolaryngol Head Neck Surg* 141(5):442–450. <https://doi.org/10.1001/jamaoto.2015.129>
- Dietz A, Wüstefeld M, Niskanen M, Löppönen H (2016) cochlear implant surgery in the elderly: the feasibility of a modified supra-aural approach under local anesthesia. *Otol Neurotol* 37(5):487–491. <https://doi.org/10.1097/MAO.0000000000001025>
- European Association of Cochlear Implant Users (2017) What is the cochlear implant? <https://eurocin.eu/what-is-the-ci->. Accessed 5 Jan 2020
- Rinck P (2018) *Magnet Resonance in Medicine. A Critical Introduction*, 12th edn. www.rinckside.org
- Statistisches Bundesamt (2018) *Fallpauschalenbezogene Krankenhausstatistik (DRG-Statistik), Diagnosen und Prozeduren der vollstationären Patientinnen und Patienten in Krankenhäusern*. Accessed 5 July 2018
- Kim B, Kim J, Park J, Kim S, Kim H, Choi J (2015) Adverse events and discomfort during magnet resonance imaging in Cochlear implant recipients. *Jama Otolaryngol Head Neck Surg* 141(1):45–52. <https://doi.org/10.1001/jamaoto.2014.2926>
- Carlson ML, Neff BA, Link MJ, Lane JJ, Watson RE, McGee KP, Bernstein MA, Driscoll CL (2015) Magnet resonance imaging with cochlear implant magnet in place: safety and imaging quality. *Otol Neurotol* 36(6):965–971. <https://doi.org/10.1097/MAO.0000000000000666>
- Azadarmaki R, Tubbs R, Chen DA, Shellock FG (2013) MRI information for commonly used otologic implants. *Otolaryngol Head Neck Surg* 150(4):512–519. <https://doi.org/10.1177/0194599813518306>
- Posner D, Scott A, Polite C, Lustig LR (2010) External magnet displacement in cochlear implants: causes and management. *Otol Neurotol* 31(1):88–93. <https://doi.org/10.1097/MAO.0b013e3181c2a0c4>
- Majdani O, Leinung M, Rau T, Akbarian A, Zimmerling M, Lenarz M, Lenarz T, Labadie R (2008) Demagnetization of cochlear implants and temperature changes in 3.0T MRI environment. *Otolaryngol Head Neck Surg* 139(6):833–839. <https://doi.org/10.1016/j.otohns.2008.07.026>
- Shew M, Wichova H, Lin J, Ledbetter LN, Staecker H (2018) Magnet resonance imaging with cochlear implants and auditory brainstem implants: are we truly practicing MRI safety? *Laryngoscope* 129(2):482–489. <https://doi.org/10.1002/lary.27516> (Epub 2018 Nov 9)
- Young NM, Rojas C, Deng J, Burrows D, Ryan M (2016) Magnetic resonance imaging of cochlear implant recipients. *Otol Neurotol* 37(6):665–671. <https://doi.org/10.1097/MAO.0000000000001053>
- Hassepass F, Stabenau V, Maier W, Arndt S, Laszig R, Beck R, Aschendorff A (2014) Revision surgery due to magnet dislocation in cochlear implant patients: an emerging complication. *Otol Neurotol* 35(1):29–34. <https://doi.org/10.1097/MAO.0b013e3182a5d2c5>
- Leong W, Yuen H (2018) Dislocation of cochlear implant magnet during 1.5 Tesla magnet resonance imaging despite head bandaging, and its repositioning using an endoscopic approach. *J Laryngol Otol* 132(10):943–945. <https://doi.org/10.1017/S0022215118001421> (Epub 2018 Aug 28)
- Shew M, Bertsch J, Camarata P, Staecker H (2017) Magnet resonance imaging in a neurofibromatosis type 2 patient with a novel MRI-compatible auditory brainstem implant. *J Neurol Surg Rep* 78(1):e12–e14. <https://doi.org/10.1055/s-0036-1597588>
- Todt I, Tittel A, Ernst A, Mittmann P, Murze S (2017) Pain free 3 T MRI scans in Cochlear implantees. *Otol Neurotol* 38(10):e401–e404
- Nospes S, Brockmann MA, Läßig A (2019) MRI in patients with auditory implants equipped with implanted magnets—an update: Overview and procedural management. *Radiologe*. 59(1):48–56. <https://doi.org/10.1007/s00117-018-0462-9>
- Erhardt JB, Führer E, Gruschke OG, Leupold J, Wapler MC, Hennig J, Stieglitz T, Korvink JG (2018) Should patients with brain implants undergo MRI? *J Neural Eng* 15(4):041002. <https://doi.org/10.1088/1741-2552/aab4e4>

25. Gubbels SP, McMenemy SO (2006) Safety Study of the Cochlear Nucleus® 24 Device with Internal Magnet in the 1.5 Tesla magnet resonance imaging scanner. *Laryngoscope* 116(6):865–871. <https://doi.org/10.1097/01.MLG.0000216807.03225.CE>
26. Broomfield SJ, Da Cruz M, Gibson WP (2013) Cochlear implants and magnetic resonance scans: a case report and review. *Cochlear Implants Int* 14(1):51–55. <https://doi.org/10.1179/1754762811Y.0000000027>
27. Cuda D, Murri A, Italice SG (2013) Focused tight dressing does not prevent cochlear implant magnet migration under 1.5 Tesla MRI. *Acta Otorhinolaryngol Ital* 33(2):133–136
28. Deneuve S, Loundon N, Leboulanger N, Rouillon I, Garabedian EN (2008) Cochlear implant magnet displacement during magnet resonance imaging. *Otol Neurotol* 29(6):789–190. <https://doi.org/10.1097/MAO.0b013e3181825695>
29. Advanced Bionics (2020) MRI Safety Information. <https://advancedbionics.com/de/de/home/professionals/mri-safety-information.html>. Accessed 5 Jan 2020
30. Cochlear Ltd. (2020) Cochlear™ Nucleus® Implantate und ihre MRT-Kompatibilität. <https://www.cochlear.com/de/startseite/unsere-produkte/cochlea-implantate/portfolio-der-nucleus-implantate/mrt-kompatibilitaet>. Accessed 5 Jan 2020
31. MED-EL (2019) MED-EL Hörimplantate. MRT-Sicherheitsstatus im Überblick. [https://www.medel.com/docs/default-source/isi-important-safety-information/languages/de-german/mri-safety-status-overview/aw36744_40_mri-safety-status-overview-\(de-german\).pdf?sfvrsn=7f126942_2](https://www.medel.com/docs/default-source/isi-important-safety-information/languages/de-german/mri-safety-status-overview/aw36744_40_mri-safety-status-overview-(de-german).pdf?sfvrsn=7f126942_2). Accessed 5 Jan 2020
32. Oticon Medical (2020) MRT und Cochlea-implantate. https://wdh01.azureedge.net/de/-/media/medical/main/files/for-professionals/ci/surgical/mri/en/210437uk_qg_mri-checklist_version_a_2019-07_low.pdf?la=en&rev=7D8E&hash=2FABD0D66C26CF4D05602040333F0CCA. Accessed 5 Jan 2020
33. Nospes S, Mann W, Keilmann A (2013) Magnetic resonance imaging in patients with magnetic hearing implants: overview and procedural management. *Radiologe*. 53(11):1026–1032. <https://doi.org/10.1007/s00117-013-2570-x>
34. Hassepass F, Stabenau V, Arndt S, Beck R, Bulla S, Grauvogel T, Aschendorff A (2014) Magnet dislocation: an increasing and serious complication following MRI in patients with cochlear implants. *Rofo*. 186(7):680–685. <https://doi.org/10.1055/s-0033-1356238> (Epub 2014 Feb 4)
35. Grupe G, Wagner J, Hofmann S, Stratmann A, Mittmann P, Ernst A, Todt I (2017) Prevalence and complications of MRI scans of cochlear implant patients: english version. *HNO* 65(Suppl 1):35–40. <https://doi.org/10.1007/s00106-016-0129-7>
36. Di Nardo W, Giannantonio S, Schinaia L, De Corso E, Pauludetti G (2013) Noninvasive management of cochlear implant's inner magnet displacement after magnetic resonance: a case report. *Laryngoscope* 123:783–786
37. Özgür A, Dursuna E, Celiker FB, Terzi S (2019) Magnet dislocation during 3 T magnetic resonance imaging in a pediatric case with cochlear implant. *Braz J Otorhinolaryngol* 85(6):799–802
38. Ray J, Proops D, Donaldson I, Fielden C, Cooper H (2004) Explantation and reimplantation of cochlear implants. *Cochlear Implants Int* 5:160–167

Publisher's Note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

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1. von Ilberg CA, Baumann U, Kiefer J, Tillein J, Adunka OF (2011) Electric-acoustic stimulation of the auditory system: a review of the first decade. *Audiol Neurotol* 16(Suppl 2):1–30. <https://doi.org/10.1159/000327765>
2. Holden LK, Finley CC, Firszt JB, Holden TA, Brenner C, Potts LG, Gotter BD, Vanderhoof SS, Mispagel K, Heydebrand G, Skinner MW (2013) Factors affecting open-set word recognition in adults with cochlear implants. *Ear Hear* 34(3):342–360. <https://doi.org/10.1097/AUD.0b013e3182741aa7>
3. Mudery JA, Francis R, McCrary H, Jacob A (2017) Older individuals meeting medicare cochlear implant candidacy criteria in noise but not in quiet: are these patients improved by surgery? *Otol Neurotol* 38(2):187–191. <https://doi.org/10.1097/MAO.0000000000001271>
4. Kraaijenga VJC, Ramakers GGJ, Smulders YE, van Zon A, Stegeman I, Smit AL, Stokroos RJ, Hendrice N, Free RH, Maat B, Frijns JHM, Briaire JJ, Mylanus EAM, Huinck WJ, Van Zanten GA, Grohman W (2017) Objective and subjective measures of simultaneous vs sequential bilateral cochlear implants in adults: a randomized clinical trial. *JAMA Otolaryngol Head Neck Surg* 143(9):881–890. <https://doi.org/10.1001/jamaoto.2017.0745>
5. Goupell MJ, Kan A, Litovsky RY (2016) Spatial attention in bilateral cochlear-implant users. *J Acoust Soc Am* 140(3):1652
6. Dorman MF, Natale SC, Butts AM, Zeitler DM, Carlson ML (2017) The sound quality of Cochlear implants: studies with single-sided deaf patients. *Otol Neurotol* 38(8):e268–e273. <https://doi.org/10.1097/MAO.0000000000001449>
7. Mosnier I, Bebear JP, Marx M, Fraysse B, Truy E, Lina-Granade G, Mondain M, Sterkers-Artières F, Bordure P, Robier A, Godey B, Meyer B, Frachet B, Poncet-Wallet C, Bouccara D, Sterkers O (2015) Improvement of cognitive function after cochlear implantation in elderly patients. *JAMA Otolaryngol Head Neck Surg* 141(5):442–450. <https://doi.org/10.1001/jamaoto.2015.129>
8. Dietz A, Wüstefeld M, Niskanen M, Löppönen H (2016) Cochlear implant surgery in the elderly: the feasibility of a modified suprameatal approach under local anesthesia. *Otol Neurotol* 37(5):487–491. <https://doi.org/10.1097/MAO.0000000000001025>
9. European Association of Cochlear Implant Users (2017) What is the cochlear implant? <https://eurociu.eu/what-is-the-ci/>. Accessed 5 Jan 2020
10. Rinck P (2018) Magnet Resonance in Medicine. A Critical Introduction. 12th edn. www.rinckside.org
11. Statistisches Bundesamt (2018) Fallpauschalenbezogene

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12. Kim B, Kim J, Park J, Kim S, Kim H, Choi J (2015) Adverse events and discomfort during magnet resonance imaging in Cochlear implant recipients. *Jama Otolaryngol Head Neck Surg* 141(1):45–52. <https://doi.org/10.1001/jamaoto.2014.2926>
13. Carlson ML, Neff BA, Link MJ, Lane JL, Watson RE, McGee KP, Bernstein MA, Driscoll CL (2015) Magnet resonance imaging with cochlear implant magnet in place: safety and imaging quality. *Otol Neurotol* 36(6):965–971. <https://doi.org/10.1097/MAO.0000000000000666>
14. Azadarmaki R, Tubbs R, Chen DA, Shellock FG (2013) MRI Information for commonly used otologic implants. *Otolaryngol Head Neck Surg* 150(4):512–519. <https://doi.org/10.1177/0194599813518306>
15. Posner D, Scott A, Polite C, Lustig LR (2010) External magnet displacement in cochlear implants: causes and management. *Otol Neurotol* 31(1):88–93. <https://doi.org/10.1097/MAO.0b013e3181c2a0c4>
16. Majdani O, Leimung M, Rau T, Akbarian A, Zimmerling M, Lenarz M, Lenarz T, Labadie R (2008) Demagnetization of cochlear implants and temperature changes in 3.0T MRI environment. *Otolaryngol Head Neck Surg* 139(6):833–839. <https://doi.org/10.1016/j.otohns.2008.07.026>
17. Shew M, Wichova H, Lin J, Ledbetter LN, Staecker H (2018) Magnet resonance imaging with cochlear implants and auditory brainstem implants: are we truly practicing MRI safety? *Laryngoscope* 129(2):482–489. <https://doi.org/10.1002/lary.27516> (Epub 2018 Nov 9)
18. Young NM, Rojas C, Deng J, Burrowes D, Ryan M (2016) Magnetic resonance imaging of cochlear implant recipients. *Otol Neurotol* 37(6):665–671. <https://doi.org/10.1097/MAO.0000000000001053>
19. Hassepass F, Stabenau V, Maier W, Arndt S, Laszig R, Beck R, Aschendorff A (2014) Revision surgery due to magnet dislocation in cochlear implant patients: an emerging complication. *Otol Neurotol* 35(1):29–34. <https://doi.org/10.1097/MAO.0b013e3182a5d2c5>
20. Leong W, Yuen H (2018) Dislocation of cochlear implant magnet during 1.5 Tesla magnet resonance imaging despite head bandaging, and its repositioning using an endoscopic approach. *J Laryngol Otol* 132(10):943–945. <https://doi.org/10.1017/S0022215118001421> (Epub 2018 Aug 28)
21. Shew M, Bertsch J, Camarata P, Staecker H (2017) Magnet resonance imaging in a neurofibromatosis type 2 patient with a novel MRI-compatible auditory brainstem implant. *J Neurol Surg Rep* 78(1):e12–e14. <https://doi.org/10.1055/s-0036-1597588>
22. Todt I, Tittel A, Ernst A, Mittmann P, Mutze S (2017) Pain free 3 T MRI scans in Cochlear implantees. *Otol Neurotol* 38(10):e401–e404
23. Nospes S, Brockmann MA, Läßig A (2019) MRI in patients with auditory implants equipped with implanted magnets—an update: Overview and procedural management. *Radiologe* 59(1):48–56. <https://doi.org/10.1007/s00117-018-0462-9>
24. Erhardt JB, Fuhrer E, Gruschke OG, Leupold J, Wapler MC, Hennig J, Stieglitz T, Korvink JG (2018) Should patients with brain implants undergo MRI? *J Neural Eng* 15(4):041002. <https://doi.org/10.1088/1741-2552/aab4e4>

25. Gubbels SP, McMenemy SO (2006) Safety Study of the Cochlear Nucleus® 24 Device with Internal Magnet in the 1.5 Tesla magnet resonance imaging scanner. *Laryngoscope* 116(6):865–871. <https://doi.org/10.1097/01.MLG.0000216807.03225.CE>
26. Broomfield SJ, Da Cruz M, Gibson WP (2013) Cochlear implants and magnetic resonance scans: a case report and review. *Cochlear Implants Int* 14(1):51–55. <https://doi.org/10.1179/1754762811Y.0000000027>
27. Cuda D, Murri A, Italica SG (2013) Focused tight dressing does not prevent cochlear implant magnet migration under 1.5 Tesla MRI. *Acta Otorhinolaryngol Ital* 33(2):133–136
28. Deneuve S, Loundon N, Leboulanger N, Rouillon I, Garabedian EN (2008) Cochlear implant magnet displacement during magnetic resonance imaging. *Otol Neurotol* 29(6):789–190. <https://doi.org/10.1097/MAO.0b013e3181825695>
29. Advanced Bionics (2020) MRI Safety Information. <https://advancedbionics.com/de/de/home/professionals/mri-safety-information.html>. Accessed 5 Jan 2020
30. Cochlear Ltd. (2020) Cochlear™ Nucleus® Implantate und ihre MRT-Kompatibilität. <https://www.cochlear.com/de/startseite/unsere-produkte/cochlea-implantate/portfolio-der-nucleus-implantate/mrt-kompatibilitaet>. Accessed 5 Jan 2020
31. MED-EL (2019) MED-EL Hörimplantate. MRT-Sicherheitsstatus im Überblick. [https://www.medel.com/docs/default-source/isi-important-safety-information/languages/de-german/mri-safety-status-overview/aw36744_40_mri-safety-status-overview-\(de-german\).pdf?sfvrsn=7f126942_2](https://www.medel.com/docs/default-source/isi-important-safety-information/languages/de-german/mri-safety-status-overview/aw36744_40_mri-safety-status-overview-(de-german).pdf?sfvrsn=7f126942_2). Accessed 5 Jan 2020
32. Oticon Medical (2020) MRT und Cochlea-Implantate. https://wdh01.azureedge.net/de/-/media/medical/main/files/for-professionals/ci/surgical/mri/en/210437uk_qg_mri-checklist_version_a_2019-07_low.pdf?la=en&rev=7D8E&hash=2FABD0D66C26CF4D05602040333F0CCA. Accessed 5 Jan 2020
33. Nospes S, Mann W, Keilmann A (2013) Magnetic resonance imaging in patients with magnetic hearing implants: overview and procedural management. *Radiologe* 53(11):1026–1032. <https://doi.org/10.1007/s00117-013-2570-x>
34. Hasepass F, Stabenau V, Arndt S, Beck R, Bulla S, Grauvogel T, Aschendorff A (2014) Magnet dislocation: an increasing and serious complication following MRI in patients with cochlear implants. *Rofo* 186(7):680–685. <https://doi.org/10.1055/s-0033-1356238> (Epub 2014 Feb 4)
35. Grupe G, Wagner J, Hofmann S, Stratmann A, Mittmann P, Ernst A, Todt I (2017) Prevalence and complications of MRI scans of cochlear implant patients: english version. *HNO* 65(Suppl 1):35–40. <https://doi.org/10.1007/s00106-016-0129-7>
36. Di Nardo W, Giannantonio S, Schinaia L, De Corso E, Pauludetti G (2013) Noninvasive management of cochlear implant's inner magnet displacement after magnetic resonance: a case report. *Laryngoscope* 123:783–786
37. Özgür A, Dursuna E, Celiker FB, Terzi S (2019) Magnet dislocation during 3 T magnetic resonance imaging in a pediatric case with cochlear implant. *Braz J Otorhinolaryngol* 85(6):799–802
38. Ray J, Proops D, Donaldson I, Fielden C, Cooper H (2004) Explantation and reimplantation of cochlear implants. *Cochlear Implants Int* 5:160–167

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